

How to use preparative HPLC - Part 1


Evaluation of loading amount

In preparative HPLC, it is quite important to investigate loading amount per one injection because required times of injection and suitable column dimension can be estimated from loading amount.

In this note, how to evaluate loading amount using analytical column (4.6 mm I.D.) is described.

(K. Kanno)

Conditions of preparative HPLC are generally optimized as the following flowchart.



<What is important to evaluate loading amount>

Main purpose of preparative HPLC is not quantitation but purification of sample. Therefore, good peak shape is not necessarily required. It does not matter, either, if injected amount exceeds detector's upper limit of linearity. Following two points are important to evaluate loading amount;

- 1) All amount of the target compound is well retained on column**
- 2) Target compound is successfully separated from other peaks**

In following pages, how to evaluate loading amount with analytical column is shown with a focus on the two points.


1) All amount of the target compound is well retained on column

- ① Prepare sample solution at as high concentration as possible.
- ② Begin with small amount, and then increase injection volume step-by-step.

In the figure shown below, injection volume was changed as follows; 10 μ L \rightarrow 50 μ L \rightarrow 100 μ L \rightarrow 250 μ L. It should be efficient to evaluate peak shape every time increasing injection volume from 2 to 10 times.

- ③ Peak shape gets worse as injection volume is increased. Injection volume at which peak shape begins to be deteriorated is the maximum amount under the conditions.

In the figure shown below, the peak began to lose its shape when 100 μ L of sample solution was injected. Peak was split into two with 250 μ L injection. Therefore, it can be said that upper limit of injection volume (amount) under the condition should be around 100 μ L (10 mg).



In case the maximum injection volume obtained with 4.6-mm I.D. column is 100 μ L as the figure shown above, the maximum injection volume of preparative column is estimated as follows. Detail of scaling up is described in LC Technical Note No.134.

Inner diameter (I.D.) of column	4.6 mm	10 mm	14 mm	20 mm	30 mm
Relative cross-sectional area	1	4.7	9.3	18.9	42.5
The maximum injection volume	100 μ L	470 μ L	930 μ L	1.89 mL	4.25 mL


2) Target compound is successfully separated from other peaks

- ① Prepare sample solution at as high concentration as possible.
- ② Begin with small amount, and then increase injection volume step-by-step.

In the figure shown below, injection volume was changed as follows; 50 μL \rightarrow 100 μL \rightarrow 200 μL \rightarrow 400 μL . It should be efficient to evaluate separation of compounds every time increasing injection volume from 2 to 10 times.

- ③ The separation gets worse as injection volume is increased. Peak resolution (R_s) is quite useful to evaluate whether sufficient separation is achieved. Target value of R_s varies depending on required purity or how to fractionate the peak, while baseline separation ($R_s > 1.5$) is ideal. R_s value of between 1.2 and 1.5 is often sufficient when fractionation is carried by detecting peak height or peak slope.

In the figure shown below, it can be said that upper limit of injection volume (total loading amount) should be around 200 μL (18 mg) in case the target value of R_s is set at 1.2.



In case the maximum injection volume obtained with 4.6-mm I.D. column is 200 μL as the figure shown above, the maximum injection volume of preparative column is estimated as follows. Detail of scaling up is described in LC Technical Note No.134.

Inner diameter (I.D.) of column	4.6 mm	10 mm	14 mm	20 mm	30 mm
Relative cross-sectional area	1	4.7	9.3	18.9	42.5
The maximum injection volume	200 μL	940 μL	1.86 mL	3.78 mL	8.5 mL

Preparative HPLC system : PLC761




- ① G-Prep AS (autosampler)
- ② CO705C (column oven)
- ③ PU714 (pump)
- ④ UV702 (UV detector)
- ⑤ RI704P (RI detector)
- ⑥ SC762 (system controller)
- ⑦ PLC-6RcA (recycle valve)
- ⑧ FC204A (fraction collector)

Inert Family

“Inertsil” High Performance LC Column

“InertCap” GC Capillary Column

“InertSep” Solid Phase Extraction for pretreatment of analysis



Contact us

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-chome, Shinjuku-ku, Tokyo, 163-1130 Japan

TEL: +81-3 (5323)6620 FAX: +81-3 (5323)6621

GL Sciences, Inc. USA

4733 Torrance Blvd. Ste 255, Torrance, CA 90503

Tel: (310)265-4424 FAX (310)265-4425

Distributors Outside of Japan and USA

GL Sciences uses distributors in many countries.

You can find a local distributor in your country in the following url.

<http://www.glsciences.com/products/contact.html>

E-MAIL: world@glsc.co.jp

